



Ion Chromatography – A Reliable, Accurate Way of Measuring Ionic Pollutants in Environmental Samples

Joachim Weiss

Dionex Corporation, International Operations, Idstein, Germany

Outline

- Conventional water analysis
- Bromate analysis by conventional IC and RFIC
- Bromate analysis with post-column derivatization
- Two-dimensional analysis of bromate with suppressed conductivity detection
- Bromate and halocetic acid analysis by IC-MS
- Summary

Advantages of Ion Chromatography

- **High Sensitivity Detection**

Allows sample to be diluted to reduce the concentration of matrix components

- **Specific Detection**

Analytes can be detected at low levels in the presence of much higher levels of matrix components

- **Analyte Specific Separations**

Column selectivities tailored to specific analyte classes



Conventional Water Analysis


IonPac AS22 Column Features

- Carbonate-based anion exchanger for inorganic anions in diverse sample matrices
- Targeted for AS4A-SC, AS12A, AS14, and AS14A inorganic anion applications
- Improved performance compared to AS4A-SC, AS12A, AS14, and AS14A
 - Improved peak shape and efficiencies
 - Improved selectivity
 - Carbonate resolved from inorganic anions
- High capacity
 - 210 $\mu\text{equiv/column}$ for 250 mm x 4 mm ID, super-macroporous
 - 6 $\mu\text{equiv/column}$ for 50 mm x 4 mm ID, microporous

IonPac AS22 Column Features / cont.

- Column selectivity optimized for 30°C
 - Ensures reproducible retention times
- Meets U.S. EPA Method 300.0 (A) and ISO performance requirements
- Compatible with ICS-2000 and ICS-3000 Systems
 - Generate carbonate/bicarbonate eluant with EGC- K_2CO_3 and EPM
- Compatible with ICS-90/1000/1500



Separation of Common Inorganic Anions and Acetate on IonPac AS22



Column: IonPac AG22, AS22, 4 mm
Eluant: 4.5 mmol/L Na₂CO₃/
1.4 mmol/L NaHCO₃
Temperature: 30°C
Flow rate: 1.2 mL/min
Inj. volume: 10 µL
Detection: Suppressed conductivity,
ASRS ULTRA II 4 mm,
AutoSuppression, recycle mode

Peaks	
1. Fluoride	5 mg/L
2. Acetate	20
3. Chloride	10
4. Nitrite	15
5. Bromide	25
6. Nitrate	25
7. Phosphate	40
8. Sulfate	30

Determination of Inorganic Anions in Municipal Drinking Water on IonPac AS22



Column: IonPac AG22, AS22 4 mm
 Eluant: 4.5 mmol/L Na₂CO₃/
 1.4 mmol/L NaHCO₃
 Temperature: 30°C
 Flow rate: 1.2 mL/min
 Inj. volume: 25 µL
 Detection: Suppressed conductivity,
 ASRS ULTRA II 4 mm,
 AutoSuppression, recycle mode

Peaks	1. Fluoride	0.84 mg/L
	2. Formate	0.03
	3. Chloride	15.59
	4. Nitrite	0.01
	5. Unknown	-
	6. Chlorate	0.18
	7. Bromide	0.02
	8. Nitrate	0.89
	9. Carbonate	-
	10. Phosphate	0.22
	11. Sulfate	20.45

Separation of Environmental Anions on IonPac AS22

Column: IonPac AG22, AS22, 4 mm

Temperature: 30°C

Eluant : 4.8 mmol/L Na₂CO₃/
1.0 mmol/L NaHCO₃


Inj. volume: 10 µL

Flow rate: 1.5 mL/min


Detection: Suppressed conductivity,
ASRS ULTRA II, 4 mm,
AutoSuppression,
Recycle mode

Peaks:

1. Fluoride	3 mg/L
2. Acetate	20
3. Formate	10
4. Bromate	20
5. Chloride	5
6. Nitrite	15
7. Chlorate	25
8. Bromide	25
9. Nitrate	25
10. Phosphate	40
11. Sulfate	30
12. Selenate	30
13. Iodide	40
14. Arsenate	30
15. Thiocyanate	40
16. Perchlorate	40
17. Thiosulfate	40
18. Chromate	40



Anion Analysis in High Chloride Wastewater with a High Capacity Anion Exchanger: Conductivity Detection



Column: IonPac AG9-HC, AS9-HC
Eluant: 9 mmol/L Na₂CO₃
Flow rate: 1 mL/min
Inj. volume: 25 μL
Detection: Suppressed conductivity,
ASRS-ULTRA, Recycle Mode


Peaks:	1. Fluoride	1.14	mg/L
	2. Acetate	-	
	3. Chloride	-	
	4. Nitrite	0.14	
	5. Bromide	0.82	
	6. Nitrate	1.47	
	7. Phosphate	21.1	
	8. Sulfate	30.3	

Sample: Synthetic wastewater
containing 2 g/L chloride

Anion Analysis in High Chloride Wastewater with a High Capacity Anion Exchanger : UV Detection

Column: IonPac AG9-HC, AS9-HC
Eluant: 9 mmol/L Na₂CO₃
Flow rate: 1 mL/min
Inj. volume: 25 μL
Detection: UV at 214 nm

Peaks: 1. Nitrite 0.18 mg/L
2. Bromide 0.79
3. Nitrate 1.58



Sample: Synthetic wastewater containing 2 g/L chloride

Maximum Cl:NO₂ Ratios Applicable to IonPac AS9-HC^a

Column / Detection Mode	Cl : NO ₂	Cl : NO ₂ -N
AS9-HC / Conductivity	2000:0.2 mg/L (10,000:1)	2000:0.068 mg/L (32,850:1)
AS9-HC / UV (214 nm)	5000:0.045 mg/L (111,000:1)	5000:0.014 mg/L (365,000:1)

^a Standard eluant conditions were used with a 25- μ L injection volume

EPA Method 300.0 (A) with RFIC: IonPac AS17


Column: IonPac AG17, AS17
Eluant: KOH (EG), 1 mmol/L for 1.5 min to 20 mmol/L in 3.5 min and from 20 mmol/L to 40 mmol/L in 2 min

Flow rate: 2 mL/min


Inj. volume: 25 μ L

Peaks:

1. Fluoride	2	mg/L
2. Chloride	3	
3. Nitrite	5	
4. Bromide	10	
5. Nitrate	10	
6. Carbonate	—	
7. Sulfate	15	
8. Phosphate	15	



Analysis of Sulfide and Cyanide




Separator: CarboPac PA1 with Guard
Eluant: 0.1 mol/L NaOH + 0.5 mol/L NaOAc + 0.5 % (v/v) ethylenediamine

Flow rate: 1 mL/min


Detection: DC Amperometry, Ag electrode (0 V)

Inj. volume: 50 μL (250 μL)

Sample preparation: Samples are acidified and distilled acc. to DIN 38 405 D13, and absorbed in 0.1-1 mol/L NaOH.



Separation of Alkali- and Alkaline-Earth Metals and Ammonium




Separator: IonPac CS12A
Eluant: 18 mmol/L MSA
Flow rate: 1 mL/min
Inj. volume: 25 μL
Detection: Suppressed conductivity,
Recycle Mode

Peaks:

1. Lithium	1 mg/L
2. Sodium	4
3. Ammonium	5
4. Potassium	10
5. Rubidium	10
6. Cesium	10
7. Magnesium	5
8. Calcium	10
9. Strontium	10
10. Barium	10


Cation Analysis in Municipal Drinking Water



Separator: IonPac CS12 / CG12
Eluant: 9 mmol/L H₂SO₄
Flow rate: 1 mL/min
Detection: Suppressed conductivity
Inj. volume: 25 μL
Sample: Tap water, acidified


Tap Water Analysis Utilizing a High Capacity Cation Exchanger

30 mmol/L MSA




Column: IonPac CS16 (5 µm)
(250 mm x 5 mm ID)
Flow rate: 1 mL/min
Temperature: 40°C
Inj. volume: 25 µL
Detection: Suppressed conductivity,
CSRS-ULTRA (4-mm),
AutoSuppression, Recycle Mode

48 mmol/L MSA



Peaks:	Calculated
1. Lithium	0.002 mg/L
2. Sodium	19.730
3. Ammonium	0.065
4. Potassium	0.987
5. Magnesium	7.210
6. Calcium	18.544

Separation of Chromium(III)/(VI) Utilizing Element-Specific Detection



Column: OmniPac PAX100 Guard
Eluant: 0.1 mol/L HNO₃ (pH 2 with NH₄OH)
Flow rate: 2 mL/min
Inj. volume: 250 μL
Detection: ICP-OES, 267.716 nm
Peaks: 1 mg/L Cr(III) and Cr(VI)



Bromate Analysis by Conventional IC

IonPac AS23 Column Features


- Carbonate-based anion exchange column for inorganic anions and oxyhalides
- Improved performance as compared to IonPac AS9-HC
 - Improved peak shape and efficiencies
 - Improved selectivity
 - Carbonate resolved from oxyhalides and inorganic anions
- High capacity
 - 320 μ equiv/column for 250 mm x 4 mm ID, super-macroporous
 - 6 μ equiv/column for 50 mm x 4 mm ID, microporous
 - Allows large loop injections
- Quantification of trace bromate in high ionic strength samples without sample pre-treatment or pre-concentration

IonPac AS23 Column Features / cont.

- Meets U.S. EPA Methods 300.0 (B) and 300.1 (A&B) and ISO performance requirements
- Compatible with ICS-2000 and ICS-3000 systems
 - Generate carbonate/bicarbonate eluant with EGC- K_2CO_3 and EPM
- Compatible with ICS-90/1000/1500


Determination of Oxyhalides and Inorganic Anions on IonPac AS23

Column: IonPac AG23, AS23, 4 mm
Eluant: 4.5 mmol/L Na₂CO₃/
0.8 mmol/L NaHCO₃
Temperature: 30°C
Flow rate: 1.0 mL/min
Inj. volume: 25 µL
Detection: Suppressed conductivity,
ASRS ULTRA II, 4 mm,
AutoSuppression, recycle mode



Peaks		
1. Fluoride	3	mg/L
2. Chlorite	10	
3. Bromate	20	
4. Chloride	6	
5. Nitrite	10	
6. Bromide	25	
7. Chlorate	25	
8. Nitrate	25	
9. Phosphate	40	
10. Sulfate	30	

Trace Analysis of Bromate with IonPac AS23/CRD and Suppressed Conductivity Detection




Column: AS23, AG23
Eluant: 4.5 mmol/L Na_2CO_3
1.0 mmol/L NaHCO_3
Flow rate: 1.0 mL/min
Detection: Suppressed conductivity (Atlas + 2 CRD)
Scrubber: 0.5 mol/L KOH (Recycle Mode)
Background: < 1.0 μS
Inj. volume: 100 μL

<u>Peak</u>	<u>ppb</u>
Fluoride	10
Chlorite	10
Bromate	5
Chloride	100
Chlorate	10
Bromide	500
Nitrate	500
Sulfate	100

IonPac AS19 Column Features

- Hydroxide-selective anion exchanger for trace analysis of bromate in drinking water
- Use of RFIC for an environmental key application
- High ion-exchange capacity
 - 240 $\mu\text{equiv/column}$ for 250 mm x 4 mm ID, super-macroporous
 - 6 $\mu\text{equiv/column}$ for 50 mm x 4 mm ID, microporous
 - Compatible with large-volume injection
- Quantitation of bromate at trace level in samples with high ionic strength without sample preparation or pre-concentration

Trace Analysis of Bromate on IonPac AS19 with RFIC




Column: Ion Pac AG19/AS19 (4-mm)
Flow rate: 0.75 mL/min
Eluant: 20 mmol/L KOH (EG)
Detection: Suppressed conductivity / External Water Mode
Inj. volume: 100uL
Sample: 2 ppb bromate in 100 ppm chloride and sulfate





Bromate Analysis with Post-column Derivatization

System Configuration for EPA Method 326.0



Low Level DBP Anion Standard Using EPA Method 326.0



Column: IonPac AG9-HC, AS9-HC
Eluant: 9 mmol/L Na₂CO₃
Temperature: 30°C
Flow rate: 1.3 mL/min
Inj. volume: 225 μL
Detection: A) Suppressed conductivity,
External water mode

B) UV (352 nm)

Post-column reagent: Acidified KI

PCR flow rate: 0.4 mL/min


Post-column heater: 80°C

Peaks:

1. Chlorite	5 μg/L
2. Bromate	5
3. DCA*	
4. Bromide	5
5. Chlorate	5

* DCA = Dichloroacetate
Quality control surrogate

Analysis of Drinking Water Using EPA Method 326.0




Column: IonPac AG9-HC, AS9-HC
 Eluant: 9 mmol/L Na₂CO₃
 Temperature: 30 °C
 Flow rate: 1.3 mL/min
 Inj. volume: 225 µL
 Detection: A) Suppressed conductivity, External Water Mode
 B) UV (352 nm)

Post-column reagent: Acidified KI
 PCR flow rate: 0.4 mL/min
 Post-column Heater: 80°C

Peaks:
 (A)
 1. Chlorite 1.85 µg/L
 2. Bromate 0.33
 3. DCA*
 4. Bromide 20.3
 5. Chlorate 82.9
 (B)
 2. Bromate 0.21 µg/L

Analysis of Ozonated Bottled Water Using EPA Method 326.0



Column: IonPac AG9-HC, AS9-HC
 Eluant: 9 mmol/L Na_2CO_3
 Temperature: 30 °C
 Flow rate: 1.3 mL/min
 Inj. volume: 225 μL
 Detection: A) Suppressed conductivity,

B) UV(352 nm)

Post-column reagent: Acidified KI
 PCR flow rate: 0.4 mL/min
 Post-column heater: 80 °C
 Peaks:

(A) 1. Chlorite	n.d.	$\mu\text{g/L}$
2. Bromate	1.52	
3. DCA*		
4. Bromide	1.12	
5. Chlorate	1.08	

(B) 2. Bromate 1.84 $\mu\text{g/L}$

* DCA = Dichloroacetate
 Quality control surrogate



Two-dimensional Analysis of Bromate with Suppressed Conductivity Detection


Requirements for Trace Analysis of Bromate in Drinking Water

- High resolution between analyte ions
 - Minimal interferences by matrix ions
 - Suitable for different concentrations of matrix ions
- Sensitive detection – preferably in the low $\mu\text{g/L}$ range – with suppressed conductivity detection
- Simple and direct injection technique
- Automation
 - Avoid labor-intensive SPE methods
 - Inline method
- Solution
 - Matrix diversion with a 2-D technique with hydroxide eluants
 - » High selectivity
 - » High sensitivity with suppressed conductivity detection

Instrumental Setup for Bromate Analysis

1. Dimension



2. Dimension



Features of the New 2-D Method

- Allows large-volume injection in the 1. dimension (4-mm column)
 - Possibility of large-volume injection in comparison to the standard method, because capacity and selectivity of the analytical column in the 1. dimension determines recovery; analyte analysis is carried out in the 2. dimension.
- Focuses the analytes in the concentrator column after suppression in the 1. dimension
 - Hydroxide eluant is converted to water, ideal for focussing or pre-concentrating analytes
- Sensitivity increase by using a 2-mm column with the lower flow rate in the 2. dimension
- Possibility of using separator columns with different selectivity in the 2. dimension
 - Higher selectivity
- Easy implementation with the ICS-3000

2-D Gradient Elution of Bromate on IonPac AS19



Column: A. IonPac AG19, AS19, 4 mm
 B. IonPac AG19, AS19, 2 mm

Eluant: KOH (EG)

Flow rate: A. 1.0 mL/min
 B. 0.25 mL/min

Suppressor: A. ASRS ULTRA II, 4 mm
 B. ASRS ULTRA II, 2 mm

Inj. volume: 500 µL

Concentrator: TAC-ULP1

Temperature: 30°C

Time	Gradient (mM)	Gradient (mM)
0	10	10
10	10	19,5
10	25	45
30	45	30,1
10	34,5	45
35	10	45

Peaks:	A	B
Bromate	0.005 mg/L	0.005 mg/L
Chloride	250	0.030
Sulfate	250	0.150

Bromate Recovery

Matrix Concentration [mg/L]	Bromate Peak Area	Recovery [%]
0	0.0248	100.0
50	0.0245	98.8
100	0.0250	100.8
150	0.0244	98.4
200	0.0249	100.4
250	0.0249	100.4

Sensitivity


Dimension	Peak Area	Flow Rate [mL/min]	Sensitivity
First (4 mm)	0.0063	1	1
Second (2 mm)	0.0248	0.25	3.936

Detection limit (MDL) for bromate, based on $n = 7$ runs (Student t Test), was calculated to be 0.2 ppb for a 500- μ L injection with 200 ppm chloride and sulfate using suppressed conductivity detection.

The background of the slide is a light blue, semi-transparent image. On the left side, there is a 200 mL Erlenmeyer flask with a white grid pattern on its surface. To the right of the flask is a white laboratory rack containing several test tubes. In the upper left corner, there is a ball-and-stick molecular model with blue and white spheres connected by grey rods. The overall scene is a laboratory setting.

Bromate and Haloacetic Acid Analysis with IC-MS



IC-MS Flow Diagram



ELMO = Enhanced Low Mass Option

MSQ is a trademark of ThermoFinnigan Corp.

Determination of Bromate in Drinking Water by IC-MS



Column: IonPac AG9-HC, AS9-HC (2-mm ID)
Eluant: 9 mmol/L Na₂CO₃
Flow rate: 0.25 mL/min
Inj. volume: 50 μL
Detection: A) Suppressed conductivity,
ASRS-ULTRA,
AutoSuppression
external water mode

B) Mass spec: -ESI mode, SIM,
m/z = 127, Source: 10 V
ESI probe: 275 °C, -2.5 kV
Boxcar smoothing, 15 points

Peaks: 1. Bromate 22 μg/L
2. DCAA NQ

DCAA = Dichloroacetic acid

U.S. EPA Regulation for Haloacetic Acids


- Current: Stage 1 DBPR
 - Monitor of HAA5 at all plants that disinfect with chlorine
 - » Total MCAA, MBAA, DCAA, DBAA, and TCAA
 - » Maximum Contamination Level (MCL) = 60 $\mu\text{g/L}$
 - Maximum Contamination Level Goal (MCLG)
 - » DCAA should not be present
 - » TCAA less than 30 $\mu\text{g/L}$
- Future
 - MCL less than 30 $\mu\text{g/L}$ for HAA5 proposed in Stage 2 DBPR

Development of a New IC-MS Method

- New IC-MS method
 - Separation of 9 HAA's and bromate in drinking water matrix
 - No derivatization or sample pretreatment
 - Fully automated
 - Low MDL and high recoveries at high concentration of drinking water matrix
- Optimize ESI-MS conditions
 - Cone voltage: different for each analyte
 - Temperature
- Selection of an internal standard
 - Stable
 - No interference with analytes
- Selection of a suitable separator column
- Detection limits and recoveries


Flow Diagram with Matrix Diversion

MD Ports 3, 4, and 5 are plugged



In this configuration the AXP-MS delivers acetonitrile to the MSQ continuously. The matrix diversion valve is used to divert sample matrix to waste and then send the analytical stream to the MSQ. The analytical stream is mixed with solvent in a high-efficiency static-mixing tee before entering the MSQ.

Hardware Setup for IC/MS




Column Evaluation for Ion Chromatography

- Columns investigated: AS16, AS19, and AS20
 - AS16: inadequate separation
 - AS19:
 - » Poor recovery and high MDL in presence of matrix
 - » Co-elution of sulfate, DBCAA, and TCAA
 - » Very narrow window for matrix diversion
 - AS20:
 - » Similar performance to AS19
 - » Faster separation
 - » Low capacity: peak broadening in matrix; high MDL for MBAA
- Prototype column:
 - » Significantly higher capacity
 - » Large enough time window for chloride diversion: allows resolution at 250 ppm chloride
 - » Separation of SO_4^{2-} , TCAA, and DBCAA => able to separate nine HAAs

Separation of Haloacetic Acids on IonPac AS24

Column: IonPac AG24, AS24 2 mm
 Eluant: EGC II KOH cartridge
 Temperature: 15 °C
 Flow rate: 0.3 mL/min
 Inj. volume: 100 µL

Suppressor: ASRS MS, 2 mm
 Detector: API 2000, MRM mode
 Solvent: MeCN, 0.20 mL/min
 KOH Gradient: Conc. (mM) Time (min)
 7.0 0.0
 7.0 18.0
 18.0 36.5
 60.0 36.8
 60.0 52.8
 7.0 53.0
 7.0 60.0



Peaks:		µg/L
1.	Chloroacetic acid	1
2.	Bromoacetic acid	1
3.	Dichloroacetic acid	1
4.	Bromochloroacetic acid	1
5.	Dibromoacetic acid	1
6.	Trichloroacetic acid	1
7.	Bromodichloroacetic acid	1
8.	Chlorodibromoacetic acid	1
9.	Tribromoacetic acid	1

MDL Calculations for HAA and Bromate in DI Water and Simulated Matrix

Analyte (m/z)	D.I. H ₂ O		*Simulated Matrix	
	MDL n=7 ng/L	%RSD	MDL n=7 ng/L	%RSD
MCAA (93)	110	6.4	192	5.3
DCAA (127)	88	6.2	207	9.8
Bromate (127)	119	6.6	282	12.4
MBAA (137)	271	14.5	352	16.3
TCAA (161)	307	15.6	397	23.2
BCAA (173)	210	9.8	259	8.5
DBAA (216.8)	177	10.1	262	9.3
DCBAA (206.8)	762	17.5	783	18.2
DBCBA (206.8)	837	18.3	844	12.2
TBAA (251)	894	21.1	952	31.5

* Simulated Matrix: Sulfate 250 ppm; Chloride 250 ppm; Nitrate 20 ppm; NH₄Cl 100 ppm

Recovery

Analyte	% Recovery			
	500 ppt	750 ppt	1000 ppt	2000 ppt
MCAA	97.3	98.2	96.8	99.9
DCAA	93.2	97.6	97.4	101.4
Bromate	101.0	97.3	98.5	105.2
MBAA	90.7	91.2	96.5	97.8
TCAA	92.6	94.7	98.1	102.2
BCAA	94.8	97.2	103.7	98.4
DBAA	93.1	95.2	97.6	103.5
DCBAA	B.D.	B.D.	93.5	107.2
DBCBA	B.D.	B.D.	109.7	97.9
TBAA	B.D.	B.D.	88.4	93.5

B.D.—Below Detection

Recovery was determined in a simulated sample containing 250 ppm chloride, 250 ppm sulfate, 20 ppm nitrate as compared to DI water.

Conclusions

- AS22 for conventional water analysis
- AS9-HC is ideal for high ionic strength samples
- AS17/19 for water analysis by RFIC
- High capacity cation exchangers for large Na/NH₄ concentration ratios
- Element-specific detection for chromium speciation
- Conventional IC and RFIC for DBPs
- IC-MS as an alternative to post-column derivatization and large volume injection for DBPs



Thank you

for your attention