

# Advanced technics in the determination of sum parameters (TIC/TOC/TNb)



Andreas Martens,  
Elementar Analysensysteme GmbH,  
Hanau, Germany

БЕОГРАД, 02.09.2008

What is a „sum parameter“?

## Definition: **sum parameter**

A composite parameter achieved by the determination of a general characteristic ( $x$ ) of various chemical compounds ( $x_i$ ) as sum without specification of a single compound.

$$\sum_{i=1}^i x_i = x_1 + x_2 + \dots + x_i$$

# sum parameters

...are commonly used in water and waste water analyses for the determination of its contamination with organic compounds.

Typical sum parameters are:

- BOD<sub>x</sub>** - **biological oxygen demand**
- COD** - **chemical oxygen demand**
- TOC** - **total organic carbon**

## sum parameters

**BOD<sub>x</sub>** and **CSB** are indirect methods to determine the amount of organic carbon by the determination of the biological and chemical oxygen demand during the decomposition of its compounds.

Disadvantage:

- **slow**
- **insufficient or wrong oxidation**
- **depending from biological activity**
- **heavy metal containing reagents**
- **no automation**


# sum parameters

## 1. Bacteria




Oxygen saturated bacteria solution is added to the sample

## 2.




Determination of the oxygen contents

## 3. T = const.



The closed sample is kept at 20°C. A part of the organic compounds are metabolized by bacteria under consumption of dissolved oxygen.

## 4.



After 5 days the oxygen concentration is determined. The BOD<sub>5</sub> in mg/L O<sub>2</sub> is calculated from the consumed oxygen.


**BOD<sub>5</sub> - biological oxygen demand after 5 days**

DIN 38 409 H51  
DIN 38 409 H52


**Disadvantage:**  
particulates are not sufficiently oxidized

Reference: Fonds der chemischen Industrie - Umweltbereich Wasser



# sum parameters



Sulfuric acid, silver ions (catalyst) and potassium dichromate (oxydation agent) are added to the sample.



The solution is refluxed at 148°C (boiling point) for 2 hours to oxidized dissolved organic compounds.



After cooling unconsumed Dichromate is determined by back titration with  $Fe^{2+}$ . The chemical oxygen demand is calculated from the difference with the consumed oxydation agent.

## COD - chemical oxygen demand

### DIN 38 409 H41

**Disadvantage:** higher or fluctuating results, if dichromate is consumed due to the oxidation of sulfur containing compounds

Reference: Fonds der chemischen Industrie - Umweltbereich Wasser

# sum parameters

## **EN 1484 - European Standard**

Water analysis -

Guidelines for the determination of the total organic carbon (TOC) and dissolved organic carbon (DOC)

**The EN 1484 approves the determination of TOC and DOC as appropriate method to determine the contamination of water with organic compounds.**



# sum parameters

## EN 1484 - European Standard

Determination of TOC and DOC by:

- **Wet chemical - UV-persulfat method**
- **Combustion method**

# Determination of TOC and DOC

## Wet chemical - UV-persulfat method

### Advantage:

- large sample volume (up to 20mL)
- high sensitivity

### Disadvantage:

- insufficient oxidation of particulates
- not applicable for electrolytes
- not for highly contaminated samples
- not for solid samples

# Determination of TOC and DOC

## Combustion method

### Advantage:

- **quantitative oxidation of all carbon containing compounds**
- **low matrix influence**
- **applicable for waste water and solids**

### Disadvantage:

- **lower injection volume (up to 3mL)**
- **problematically with salts**
- **lower sensitivity**

# Conclusion

The **combustion method** has the largest potential to be state of the art, when the disadvantages can be minimized.

## EN 1484 - European Standard

**TC** - total carbon

**TIC** - total inorganic carbon („H<sub>2</sub>CO<sub>3</sub>“)

**TOC** - total organic carbon (+ CN<sup>-</sup>, OCN<sup>-</sup>, SCN<sup>-</sup>)

**TIC + TOC = TC** (additive determination)

**TOC = TC - TIC** (difference method)

Direct method:

**TOC = (TC - TIC)** after external acidification

## EN 1484 - European Standard

TOC - total organic carbon

NPOC - non-purgable organic carbon

POC (VOC) - purgable organic carbon

$$\text{TOC} = \text{POC} + \text{NPOC}$$

$$\Rightarrow \text{TIC} + \text{POC} + \text{NPOC} = \text{TC}$$

or

$$\Rightarrow \text{NPOC} = \text{TC} - \text{TIC} - \text{POC}$$

$$\text{POC} = 0 \quad \Rightarrow \quad \text{NPOC} = \text{TOC}$$


# Additive determination



$$\text{TIC} + \text{TOC} = \text{TC}$$


or

$$\text{TIC} + \text{POC} + \text{NPOC} = \text{TC}$$




liqui TOC II  
liqui TOC trace

# NPOC direct (after external acidification)







# TIC/NPOC - determination



# TIC/POC/ NPOC - determination



# Dynamic 2-zone-furnace



Determination of  
TIC and TOC  
in solids  
**without  
acidification**  
by the  
**temperature  
ramp method.**

# Difference method

$$\text{TOC} = \text{TC} - \text{TIC}$$

or


$$\text{NPOC} = \text{TC} - \text{TIC} - \text{POC}$$



vario TOC cube

# vario TOC cube

## Carrier Gas Flow



Funktionssymbole  
functional icons

- Zweiwegeventil  
two way valve
- Mehrwegeventil  
multi way valve
- Rückschlagventil  
check valve
- Dosierpumpe  
dosing pump
- Flussregler  
mass flow controller

- Flusssensor  
flow sensor
- Drucksensor  
pressure sensor
- Filter  
filter
- Abfall  
waste
- Trocknung  
gas drying

- Verbrennung  
combustion
- Messzelle  
detector cell
- Magnetrührer  
magnetic stirrer
- Halogenabsorber  
halogen absorber


Basis Symbole  
basic icons

- Sensor  
sensor
- Aktuator  
actuator
- Schlauch  
plastic tube

## TIC/TC analysis

# vario TOC cube

## Acid Dosing



Funktionssymbole  
functional icons

- Zweiwegeventil  
two way valve
- Mehrwegeventil  
multi way valve
- Rückschlagventil  
check valve
- Dosierpumpe  
dosing pump
- Flussregler  
mass flow controller

- Flusssensor  
flow sensor
- Drucksensor  
pressure sensor
- Filter  
filter
- Abfall  
waste
- Trocknung  
gas drying

- Verbrennung  
combustion
- Messzelle  
detector cell
- Magnetrührer  
magnetic stirrer
- Halogenabsorber  
halogen absorber


Basis Symbole  
basic icons

- Sensor  
sensor
- Aktuator  
actuator
- Schlauch  
plastic tube

## TIC/TC analysis

# vario TOC cube

Flushing with sample



Funktionssymbole  
functional icons

- Zweiwegeventil  
two way valve
- Mehrwegeventil  
multi way valve
- Rückschlagventil  
check valve
- Dosierpumpe  
dosing pump
- Flussregler  
mass flow controller

- Flusssensor  
flow sensor
- Drucksensor  
pressure sensor
- Filter  
filter
- Abfall  
waste
- Trocknung  
gas drying

- Verbrennung  
combustion
- Messzelle  
detector cell
- Magnetrührer  
magnetic stirrer
- Halogenabsorber  
halogen absorber


Basis Symbole  
basic icons

- Sensor  
sensor
- Aktuator  
actuator
- Schlauch  
plastic tube

## TIC/TC analysis





# vario TOC cube





## Dosing (TIC-Measurement)






Funktionssymbole  
functional icons

-  Zweiwegeventil  
two way valve
-  Mehrwegeventil  
multi way valve
-  Rückschlagventil  
check valve
-  Dosierpumpe  
dosing pump
-  Flussregler  
mass flow controller

-  Flusssensor  
flow sensor
-  Drucksensor  
pressure sensor
-  Filter  
filter
-  Abfall  
waste
-  Trocknung  
gas drying

-  Verbrennung  
combustion
-  Messzelle  
detector cell
-  Magnetrührer  
magnetic stirrer
-  Halogenabsorber  
halogen absorber

Basis Symbole  
basic icons


-  Sensor  
sensor
-  Aktuator  
actuator
-  Schlauch  
plastic tube

## TIC/TC analysis



# vario TOC cube

## TIC-Measurement



Funktionssymbole  
functional icons



Zweiwegeventil  
two way valve



Mehrwegeventil  
multi way valve



Rückschlagventil  
check valve



Dosierpumpe  
dosing pump



Flussregler  
mass flow controller



Flusssensor  
flow sensor



Drucksensor  
pressure sensor



Filter  
filter



Abfall  
waste



Trocknung  
gas drying



Verbrennung  
combustion



Messzelle  
detector cell



Magnetrührer  
magnetic stirrer



Halogenabsorber  
halogen absorber

Basis Symbole  
basic icons



Sensor  
sensor



Aktuator  
actuator




Schlauch  
plastic tube

TIC/TC analysis

# vario TOC cube

## TC-Measurement



Funktionssymbole  
functional icons



Zweiwegeventil  
two way valve



Mehrwegeventil  
multi way valve



Rückschlagventil  
check valve



Dosierpumpe  
dosing pump



Flussregler  
mass flow controller



Flusssensor  
flow sensor



Drucksensor  
pressure sensor



Filter  
filter



Abfall  
waste



Trocknung  
gas drying



Verbrennung  
combustion



Messzelle  
detector cell



Magnetrührer  
magnetic stirrer



Halogenabsorber  
halogen absorber

Basis Symbole  
basic icons



Sensor  
sensor



Aktuator  
actuator




Schlauch  
plastic tube

TIC/TC analysis

# vario TOC cube

TIC/TC-Waste



Funktionssymbole  
functional icons



Zweiwegeventil  
two way valve



Mehrwegeventil  
multi way valve



Rückschlagventil  
check valve



Dosierpumpe  
dosing pump



Flussregler  
mass flow controller



Flusssensor  
flow sensor



Drucksensor  
pressure sensor



Filter  
filter



Abfall  
waste



Trocknung  
gas drying



Verbrennung  
combustion



Messzelle  
detector cell



Magnetrührer  
magnetic stirrer



Halogenabsorber  
halogen absorber

Basis Symbole  
basic icons



Sensor  
sensor



Aktuator  
actuator




Schlauch  
plastic tube

TIC/TC analysis

# vario TOC cube

Flushing with sample



Funktionssymbole  
functional icons

- Zweiwegeventil  
two way valve
- Mehrwegeventil  
multi way valve
- Rückschlagventil  
check valve
- Dosierpumpe  
dosing pump
- Flussregler  
mass flow controller

- Flusssensor  
flow sensor
- Drucksensor  
pressure sensor
- Filter  
filter
- Abfall  
waste
- Trocknung  
gas drying

- Verbrennung  
combustion
- Messzelle  
detector cell
- Magnetrührer  
magnetic stirrer
- Halogenabsorber  
halogen absorber

Basis Symbole  
basic icons

- Sensor  
sensor
- Aktuator  
actuator
- Schlauch  
plastic tube

NPOC analysis

# vario TOC cube

combines for the very first time

**110 years of experience in elemental analysis**  
and more than  
**30 years of experience in TOC analysis**



# vario TOC cube

THE POWER OF EXPERIENCE



Sets new standards  
in TOC analysis

**measuring range for  
TOC from  
2 ppb to 60.000 ppm**



**optional up to 50mg C<sub>abs.</sub> in solids**

**Makes „auto-dilution“ unnecessary.**



# vario TOC cube

THE POWER OF EXPERIENCE



Sets new standards  
in TOC analysis


with new detectors –  
free of spectral  
interferences, e.g. of  $\text{SO}_2$   
and with extraordinary  
linearity over the entire  
detection range.

# vario TOC cube

THE POWER OF EXPERIENCE

Sets new standards  
in TOC analysis

New NDIR detectors with  
outstanding  
performance in  $TN_b$  –  
determination.



Also available for  $TN_b$  determination:

**CLD 0.2ppm – 200ppm**  
**ECD 0.2ppm – 500ppm**

Working range:  
**0.05ppm – 50,000ppm**



# vario TOC cube

THE POWER OF EXPERIENCE



Sets new standards  
in TOC analysis

**with its integrated  
automatic sample  
feeder for liquids**

# vario TOC cube

THE POWER OF EXPERIENCE



Sets new standards in  
TOC analysis

**with its revolutionary  
automated analysis of  
solids**

available with

- 120 positions for up to 300 mg
- 80 positions for up to 1 g
- also applicable for liquids and suspensions by capsule sealing technique

# vario TOC cube

THE POWER OF EXPERIENCE




Sets new standards  
in TOC analysis

**with its unique reusable  
ash/salt crucible for  
complete matrix separation,  
e.g. no restrictions on salt  
content in liquids and solids**

# vario TOC cube

liquid samples



solid samples

Combustion tubes

# vario TOC cube

Combustion tube for solids



Combustion tube for liquids

Combustion tubes

# vario TOC cube

THE POWER OF EXPERIENCE



Sets new standards  
in TOC analysis

**with no restrictions  
on the particle size  
in liquids**

No need any more for inaccurate  
suspension method

# vario TOC cube

THE POWER OF EXPERIENCE

Applications from liquids to...



- sea water
- saline solutions
- electrolytes
- ...



- reagent water
- drinking water
- waste water
- ...



# vario TOC cube

THE POWER OF EXPERIENCE

... to solids

- particles
- suspensions
- slurries ...



- soils
- residues
- solids ...



# vario TOC cube


THE POWER OF EXPERIENCE

Approved methods, once again improved with IDS-3\*



condensation

counter  
gas  
membrane  
drying




chemical  
drying



\* 3-step intensive drying system

# vario TOC cube

THE POWER OF EXPERIENCE




Comprehensive  
math data  
processing



- **Statistics**
- **Manual and automatic peak integration**
- **Base line correction**
- **Peak subtraction**

# vario TOC cube

THE POWER OF EXPERIENCE




No.	Hole Pos.	Name	TIC Area	TIC [mg/L]	NPOC Area	Dil. Factor	NPOC [mg/L]
1	1	Std (1mL)	0.00	0.00	0.20	1.00	5.00
2	2	Std (1mL)	0.00	0.00	0.40	1.00	10.00
3	3	Std (1mL)	0.00	0.00	0.60	1.00	15.00
4	4	Std (1mL)	0.00	0.00	0.80	1.00	20.00
5	5	Std (1mL)	0.00	0.00	1.00	1.00	25.00

Automatic Calibration  
by fixed volume  
and different conc.

Dilution series




Automatic outlier elimination



# vario TOC cube


THE POWER OF EXPERIENCE

Automatic Calibration  
by fixed concentration  
and different volumes  
auto-dilution-effect  
without dilution errors



NPOC calibration, 25 ppm

Sample Pane (Text View)						
No.	Hole Pos.	Name	TIC Area	TIC [mg/L]	NPOC Area	Dil. Factor
1	1	Std{0.20}	0.00	0.00	0.20	1.00
2	1	Std{0.40}	0.00	0.00	0.40	1.00
3	1	Std{0.60}	0.00	0.00	0.60	1.00
4	2	Std{0.80}	0.00	0.00	0.80	1.00
5	2	Std{1.00}	0.00	0.00	1.00	1.00



Zoom factor: 1  
a = -4.097819e-010  
b = +2.9500000e-002  
r = 1.000000  
sy = 0.000000



Multi-point-calibration from a single standard solution



# vario TOC cube

THE POWER OF EXPERIENCE

Analysis of well water for animals  
(liquid samples, particles)




Method	Sample	TC [mg/L]
mean [%]	well water	76,520
RSD (absolute)		0,271
RSD (relative) %		0,354

# vario TOC cube

## THE POWER OF EXPERIENCE

Analysis of hard digestible  
silicon carbide compounds  
(SiBCN)



Method	Sample	TOC [mg/L]
mean [%]	SiBCN	30,82
RSD (absolute)		0,0683
RSD (relative) %		0,0261

vario TOC cube shows excellent oxidation behaviour even in case of hard digestible compounds like SiBCN ceramics with decomposition temperatures of approx. 2000°C.

# vario TOC cube

## THE POWER OF EXPERIENCE

Analysis of sea water  
(liquid samples, salt)



Method	Sample	TC [mg/L]
--------	--------	--------------

TC measurement of n = 50 sea water samples

mean [%]	North Sea	24,282
----------	-----------	--------

RSD (absolute)		0,160
----------------	--	-------

Method	Sample	TOC [mg/L]
--------	--------	---------------

TOC measurement of n = 50 sea water samples. External sample acidification was performed under stirring in the vials of the automatic sample feeder, followed by automatic stripping of CO<sub>2</sub> with synthetic air.


mean [%]	North Sea	12,209
----------	-----------	--------

RSD (absolute)		0,074
----------------	--	-------

**Excellent for the analysis of brine solutions**

# vario TOC cube

## THE POWER OF EXPERIENCE



### Analysis of river sediments (solid sample)

#### Method

Analysis of 25 dried river sediments  
(Elbe/Germany)  
measured in solid mode

	n	TC [%]
mean [%]	25	5,396
RSD (absolute)		0,088
Reference value (varioMAX)		5,35 %

#### Method

Analysis of 25 dried river sediments  
(Hunte/Germany)  
measured in solid mode

	n	TC [%]
mean [%]	25	4,202
RSD (absolute)		0,070
Reference value (varioMAX)		4,22 %



# vario TOC cube

THE POWER OF EXPERIENCE

## THE FACTS I

**Analytical method:** Catalytic high temperature digestion

**Standards:** ISO 8245, EPA 415.1, European standard acc. to EN 1484, ENV 12260, European Pharmacopoeia 6.0, vol.1 (20244) 2008, U.S.Pharmacopoeia-NF, USP 30 <643> 2007; DIN 38409

**Measured parameters:** TC; TOC, TIC/NPOC; DOC; POC  
optional TN<sub>b</sub>

**Digestion temperature:** free selectable up to 1200 °C,  
standard operating temperature in dependency on the  
chosen catalyst between 800 – 950 °C



# vario TOC cube

THE POWER OF EXPERIENCE

## THE FACTS II

**Measuring range (ppm):** C: 0 – 60,000 ppm (cube)

alternatively high%C

C: 0 – 100 % (up to 50mg C<sub>abs.</sub>)

C: 0 – 25,000 ppm (trace)

N: 0 – 200 ppm (CLD); 500 ppm (EC)

N: 0 – 50,000 ppm (NDIR)

**Detection limit:**

C: ± 6 µg/L SD (cube)

C: ± 3 µg/L SD (trace)

N: ± 0.02 mg/L SD (NDIR)

**Precision:**

< 1 % at > 5 mg/L C



# vario TOC cube

THE POWER OF EXPERIENCE

## THE FACTS III

**Injection volume (liquid-mode):** 0.05 – 2 mL (cube)  
0.2 – 4 mL (trace)

**Sample weight (solid-mode):**

organic substance 0.02 – 10 mg or up to 1 g soil

(depends on the sample matrix)

**Particle size:** **liquid mode:** inner tube diameter 0.8 mm

**solid mode:** free selectable in  
dependency on the chosen capsule size

**Duration of analysis:** 2 - 3 min. for TC (solids)  
3 - 4 min. per parameter (liquids)



# vario TOC cube

THE POWER OF EXPERIENCE

## THE FACTS IV

**Gas:** synthetic air (hydrocarbon  $\leq 0.1 \text{ mg}\cdot\text{L}^{-1}$   
and  $\text{CO}_2 \leq 1 \text{ mg}\cdot\text{L}^{-1}$ )

or oxygen (grade  $\geq 4.5$ )

**Gas flow:**  $200 \text{ mL}\cdot\text{min}^{-1}$  at 1 bar pre-pressure

**Dimensions:** 55 x 42 x 55 cm (L x W x H) – basic unit  
55 x 42 x 72 cm (L x W x H)

– basic unit with autosampler for liquids

**Weight:** approx. 60 kg

**Power supply:** 100/110/200/230 V, 50/60 Hz, 1.6 kW



# vario TOC cube

**THE POWER OF EXPERIENCE**

## Instrument control and data processing

Operation and control under Windows® XP or Windows® Vista. All instrument functions are digitally controllable, the comprehensive operation software includes e.g. automatic leak test, extensive fault diagnostics, monitoring of maintenance cycles, sleep- / awake function, statistical evaluation and almost unlimited memory capacity for storage of analytical data and graphics. Integration to data network and LIMS.

The possibility of remote control and diagnosis via the internet.

In full compliance with 21 CFR Part 11 (option).



# **ELEMENTAR** Analysensysteme GmbH



**Thank you for your attention.**

**ХВАЛА!**



