



# Skalar



## Analyses of inorganic parameters by Flow Analyses

Water Workshop  
University of Novi Sad  
1 – 5 september 2008



# Program

- Introduction of Skalar Analytical
- Skalar Continuous Flow Analyses
- Practical applications on CFA
  - Nutrients
  - Cyanide
  - Phenol
  - Total Nitrogen
  - Total Phosphate
  - Dissolved Organic Carbon




# Skalar Worldwide



Germany




The Netherlands



France




Belgium




**Represented in  
more than 60 countries**



Austria



United Kingdom



United States



India



Serbia



Czech Republic

# Thousands of Skalar Analyzers operational worldwide

- ✓ Commercial laboratories
- ✓ Industrial laboratories
- ✓ Universities and Research laboratories
- ✓ Process control


## Applications:

- Water
- Soil, Plant & Fertilizer
- Beer & Malt
- Food & Beverage
- Wine
- Tobacco
- Pharmaceutical
- Detergents
- Mining & Metallurgical
- Petrochemical



## According International standards

Skalar applications are according to international directives such as



- ISO
- EPA
- ASTM
- CEN
- AFNOR
- AOAC
- DIN
- NEN
- ASBC
- EBC
- Mebak
- User's methods

# Skalar Product Lines



SAN++ Analyzer  
Wet Chemistry automation



Formacs<sup>SERIES</sup> and Primacs<sup>SERIES</sup>  
TOC & TN Analyzers




Robotic Analyzers  
BOD, COD, Titrations, ISE etc



Fluo Imager  
Oil in water, Chlorofyll



Toxtracer  
Bio-assay for Toxicity



# Program


- Introduction of Skalar Analytical
- Skalar Continuous Flow Analyses
- Practical applications on CFA
  - Nutrients
  - Cyanide
  - Phenol
  - Total Nitrogen
  - Total Phosphate
  - Dissolved Organic Carbon



**SKALAR**

# Skalar Products Continuous Flow Analyzer

*San<sup>++</sup>*





# Continuous Flow Analyzer

## Typical Automated Environmental Applications

**Ammonia**

Aluminum

Boron

Bromide

Calcium

**Chloride**

Chemical Oxygen Demand

Chlorine

Cholinesterase Inhibition

Chromium VI

Color

Conductivity

**Cyanide (Total - Free - WAD) Dissolved**

**Organic Carbon**

**Fluoride**

Formaldehyde

Iron (Total - Free - Hydrolysable)

Magnesium

Manganese

Methylene Blue Active Substances (MBAS)

Molybdenum

**Nitrate + Nitrite**

Permanganate value (COD) Potassium

Total Amino Acids

Total Alkalinity

Total Carbonates

Total Hardness

**Total Nitrogen (UV and TKN)**

**Total Phosphate (UV)**

Total Phenols

Silicate

Sodium

**Sulfate**

Sulfur Dioxide

Urea


Volatile Acids

Zinc

# Continuous Flow Analyzer

## Principle of the San<sup>++</sup>

Full Instrument Control




From 1 to 16 channels  
simultaneously

Data Handling  
and Report  
Generation



# Continuous Flow Analyzer

## Modular Design



# Continuous Flow Analyzer

## Typical Advantages CFA compared to other techniques

- Implementation of complex analysis techniques like:  
inline distillation, digestion, extraction, dialysis
- No blocking by dirty samples as of mini-bore coils
- Sub ppb level sensitivity (and expended range in ppm level)
- Full end point color development (optimal chemistry)
- No synchronization required to analyze data
- No degassing of reagents

( CFA = Continuous Flow Analysis )



# Program

- Introduction of Skalar Analytical
- Skalar Continuous Flow Analyses
- Practical applications on CFA
  - Nutrients
  - Cyanide
  - Phenol
  - Total Nitrogen
  - Total Phosphate
  - Dissolved Organic Carbon



# Practical Applications on CFA

## Skalar SAN<sup>++</sup> Example Configuration

**SYSTEM 1, CONFIGURATION FOR THE SKALAR AUTOMATIC ANALYSER FOR THE ANALYSIS OF AMMONIA, SILICATE, ORTHO PHOSPHATE, NITRATE+NITRITE AND NITRITE IN WATER/SEAWATER SAMPLES**



**SA3074, sampler with diluter and mixing dials**



**SA3000, module holder**




**Lab PC with system controller, computer and printer (printer not shown)**

# Practical Applications on CFA

## Nitrite: Colorization method

According ISO 13395

Sulfanylamide



Diazonium compound



Red-purple color measured at 540 nm


**Identical to the manual method**

# Practical Applications on CFA

## Nitrite: Colorization method

Photometric detection

FLOW DIAGRAM



# Practical Applications on CFA

## Nitrate: Colorization method ( $\text{NO}_2 + \text{NO}_3$ )

According ISO 13395

$\text{NO}_3 \Rightarrow \text{NO}_2$  by reduction

Reduction methods:

- ✓ Cd-reduction
- ✓ reduction with hydraziniumsulfate
- ✓ enzymatic reduction

Sulfanylamide

↓ +  $\text{NO}_2 + \text{H}^+$

Ziazonium compound

↓  $\alpha$ -naphthylethylenediamine dihydrochloride


Red- purple color measured at 540 nm

After reduction of  $\text{NO}_3$  to  $\text{NO}_2$ :  
Identical  
measurement as  
 $\text{NO}_2$ -application

# Practical Applications on CFA

## Nitrate: Colorization method ( $\text{NO}_2 + \text{NO}_3$ )

**FLOW DIAGRAM**





# Practical Applications on CFA

## Total Cyanide, a complete Analysis Process

### Different methods:

- 1. Photometric detection with UV-A destruction**  
Total determination of CN, including thiocyanates  
According ISO 14403
- 2. Photometric detection with UV-B destruction**  
Total determination of CN, excluding thiocyanates  
According ISO 14403
- 3. Amperometric detection**  
with a silver working electrode and silver/silver chloride reference electrode  
According ISO 14403
- 4. Photometric detection of WAD (Weak and Dissociable cyanides)**  
According ISO 14403

# Practical Applications on CFA

## Total Cyanide: Photometric detection with UV-B destruction

Total determination of CN, excluding thiocyanates, according ISO 14403

CN-complex

↓ + UV-light ( + pH 3.8

HCN

↓ Separation by distillation at 125°C under vacuum

↓ Chloramine-T,

Cyanogene chloride

↓ Colourisation with 4-pyridine carboxylic acid

↓ and 1,3-dimethylbarbituric acid

Red colour, measured at 600 nm


# Practical Applications on CFA

## Total Cyanide: Photometric detection with UV-B destruction


FLOW DIAGRAM

Photometric detection

Color Reaction



In-line distillation



In-line UV digestion



Sampling

\* distillation unit is situated at the back panel of the module

\*\* glass tubing SA 5359 (2.0mm I.D., 4.0mm O.D.)

# air inlet via 50 cm 5133 tube with restrictor (3 cm 5142 tube)

## air restrictor (2CA14031)

## Cyanide, Total & Free: Method Amperometric

### **PRINCIPLE (total cyanide)**

The automated procedure for the determination of Total Cyanide is based on the following reaction: Complex bound cyanide is decomposed by the radiation of UV light, in a continuously flowing stream at pH 3.8. A UV-B lamp ( $\pm 312$  nm) and a coil of borosilicate glass is used to avoid UV light, with a wavelength of less than 290 nm, and thus preventing the conversion of thiocyanate into cyanide. The Cyanide, decomposed at pH 3.8, is kept in solution by addition of a liquid to increase the pH to  $> 10$  and then pumped to the dialyser where the sample is acidified to form hydrogen cyanide gas.

**The hydrogen cyanide gas diffuses through the hydrophobic Teflon membrane** into the alkaline receptor stream.

**The cyanide is measured amperometrically** with a silver working electrode and silver/silver chloride reference electrode. (the method is according ISO 14403)

# Practical Applications on CFA

## Phenol Index: a complete Analysis Process

Total determination of phenols, according ISO 14402

phenol-complex

↓ distillation at pH 1.4

volatile phenolic compounds

↓ potassium hexacyanoferrate(III)

Quinones


↓ 4-aminoantipyrine

yellow condensation products, measured at 505 nm

# Practical Applications on CFA

## Phenol Index: a complete Analysis Process

### FLOW DIAGRAM



# Practical Applications on CFA

## Total Nitrogen

### Different methods:

Element	Legislation	Description	Skalar system	Skalar Method
Total Nitrogen	ISO/CD 29441	Online measurement after UV-destruction	San <sup>++</sup> with online UV-destruction	M475-424
Kjeldahl-N total phosphate		Offline destruction, followed by colorimetric detection of ammonium and phosphate on CFA	San <sup>++</sup> with prior destruction on SA5640	M155-056 M503-004
Kjeldahl-N		Offline destruction, followed by colorimetric detection of ammonium on CFA	San <sup>++</sup> with prior destruction on SA5640	M155-056
Total Nitrogen	EN 12260:2003 or EN-ISO 11905-2:1997	High Temperature combustion	Formacs <sup>TM</sup>	
Totaal Nitrogen	EN-ISO 11905-1:1998	Offline destruction, followed by colorimetric detection of nitrate on CFA	San <sup>++</sup> with prior destruction on SA5640	M461-032

# Practical Applications on CFA

## Total Nitrogen by in-line UV-destruction

According ISO/CD 29441

**N-compounds**

↓ Oxidation reagent + UV-destruction

**Organic carbon to Nitrate**

↓ +T ( C )

**Ammonium to nitrate**

↓ +Cd-reduction

Identical to nitrate measurement

**Nitrate to nitrite**

↓ +sulfanilamide ,  
N-1-nafthyl-ethylen-diamine-dihydrochloride

**Red colored complex, measured at 540 nm**


Alternative for Kjeldahl



# Practical Applications on CFA

## Total Nitrogen by in-line UV-destruction

1st Order ISO 8466 Calibration - Tot.Stikstof



a = -2844.29063 b(Slope) = 35288.60938 RSD = 5760.456 r = 0.99975

# Practical Applications on CFA

## Total Nitrogen: Result comparison versus manual method

<u>Sample Identity</u>	<u>Result Skalar (TN)</u>	<u>Result ref (TN)</u>
10443292	7.5	7.5
22034434	32.9	33.7
10441790	292.2	294.0
22036038	59.8	60.6
22036040	78.8	79.3
22036043	53.2	53.1
22036149	50.5	51.8
88050901	20.9	20.6
88051101	33.0	33.4
88053301	28.9	30.3
88053501	61.1	61.9
88053601	49.3	50.6
88053602	3.8	3.7
88053701	56.8	58.0
88053801	113.4	111.0

# Practical Applications on CFA

## Total Phosphate by in-line UV-destruction

Element	Norm	Omschrijving	Skalar Systeem	Skalar Methode
Phosphate ortho & total	EN-ISO 15681-2:2003	Off-line destruction or in-line destruction	San <sup>++</sup> with in-line UV-destruction or off-line destruction	M503-004
Kjeldahl-N total phosphate		Offline destruction, followed by colorimetric detection of ammonium and phosphate on CFA	San <sup>++</sup> with prior destruction on SA5640	M155-056 M503-004

# Practical Applications on CFA

## Total Phosphate by in-line UV-destruction

According ISO 15681-2

**P-compounds**

↓ Oxidation reagent + UV-destruction

**Organic phosphate to ortho-phosphate**

↓ +T (°C)

**Inorganic phosphate to ortho-phosphate**

↓ +katalyst  
ammonium haptamolybdate


**phospho-molybdic acid complex, measured at 540 nm**

Identical to  
o-phosphate  
measurement


# Practical Applications on CFA

## Total Phosphate by in-line UV-destruction

FLOW DIAGRAM




\*\* polyethylene tube catnr. SA 5141  
\*\*\* silicone tube catnr. SA 5151



# Practical Applications on CFA

## Total Phosphate by in-line UV-destruction





# Practical Applications on CFA

## Total Phosphate / Total Nitrogen

Graph Results (UV-Digestion vv Manual)

Total Phosphate UV-P:  
Range: 1 - 500  $\mu\text{g P/l}$ .

Total Nitrogen UV-N:  
Range 20 - 5000  $\mu\text{g N/l}$ .



According to international standard regulations (ISO, DIN, etc.) accredited!


# Practical Applications on CFA

## Total Organic Carbon

### How to measure TOC?

- UV-promoted persulfate oxidation (TOC)

↳ On CFA using A TOC module, followed by IR-detection




- High temperature catalytic combustion (TOC & TN)

↳ On Formacs<sup>SERIES</sup> TOC Analyzer



# Practical Applications on CFA

## What is Total Organic Carbon?



## Practical Applications on CFA

### Total Organic Carbon – UV promoted on CFA

Sample

↓ acidified + sparged ( $N_2$ )

liberates and disperses any inorganic or volatile organic carbon  
Organic carbon rests in sample

↓ persulfate / tetra borate reagent  
+ UV digestion coil

All organic carbon oxidised to  $CO_2$


↓ acidification and sparging  
to release  $CO_2$  from solution

IR-measurement of  $CO_2$

# Practical Applications on CFA

## Total Organic Carbon – UV promoted on CFA

### FLOW DIAGRAM



# TOC / TN Analyses

## Formacs<sup>HT</sup> TOC Analyzer

New model High Temperature Combustion Analyzer



Optional  
Sampler


TOC Analyzer

Data acquisition &  
instrument control

# TOC / TN Analyses


## Formacs<sup>HT</sup>: Schematics of TOC Analyzer

### New model High Temperature Combustion Analyzer



# TOC / TN Analyses

## Formacs<sup>TN</sup>: Schematics of ND20 TN Detector



# TOC / TN Analyses

**Formacs<sup>TN</sup>: With TKN reactor**

Skalar offers a true Kjeldahl Alternative

$$TKN = TN - NN$$



- NN application offers fast and economical alternative for classical Total Kjeldahl Nitrogen analysis
- NN Reactor can also be used for IC analysis
- Low maintenance

# TOC / TN Analyses

Formacs<sup>TN</sup>: With TKN reactor

## Case Study

Italian National Research Council

Section of Hydrobiology and Ecology of Inland Waters, Verbania - Italy



For many years the water quality of the Italian lakes has been monitored on nutrients to obtain long-term trends for evaluation.

Nitrogen concentrations is an important parameter for the ecosystem functioning and needs to be monitored accurately.

Typical specifications

Range : 0.05 – 5.00 mg/l

N

CV : < 1 % F.S.

MDL : 10 µg/l N

NaNO3 Std 2.00 mg/l	TN	Results 2.00	Recovery
NaNO3 Std 4.00 mg/l	TN	4.00	%
KNO3 2 mg/l N	TN	2.04	102
NaNO2 2mg/l N	TN	2.08	104
NH4Cl 2 mg/l N	TN	2.02	101
NH4Cl 4 mg/l N	TN	4.2	105
NH4Acetate 2mg/l N	TN	2.1	105
Glycine 2 mg/l N	TN	2.07	104
Creatinine 2 mg/l N	TN	2.11	106
Creatinine 4 mg/l N	TN	4.27	107

# TOC / TN Analyses

## Formacs<sup>TN</sup> Vs Traditional TKN



- Much faster than the Kjeldahl method ; 3 minutes per measurement to 1-2 hours for TKN
- Combustion technique handles more difficult matrices (Particulate & Brines)
- No toxic chemicals.
- Many samples are measured automatically
- Better precision
- Wide dynamic range
- Elimination of operator error
  
- Conclusion: Cost reduction for working hours, chemical waste disposal and laboratory space.



**Thank you for your attention.**

**We hope to see you at the Hands-on demonstration.**

**SKALAR,**

***Your Partner in Chemistry Automation***