

Solid Phase Extraction for Water Analysis



presented by

Dr. Bernd Thierfelder

strata™ X

phenomenex®
...breaking with tradition™

Outline

1. Introduction in SPE
2. Retention Mechanisms / Sorbent Selectivity
3. strata™-X for Water Analysis
4. Applications of Water Analysis

strata™X

phenomenex®
...breaking with tradition™

1. Introduction in SPE


strata[™]X

phenomenex[®]
...breaking with tradition[™]

What is Solid Phase Extraction (SPE)?

It's a very selective and effective sample preparation method.

A day in the life of a scientist ...




Homogenization
Dilution
Weighing
Settling and Decanting
Centrifugation
Filtration
Enzymatic Hydrolysis
Evaporation
Protein Crash / Precipitation
Liquid / Liquid Extraction
Solid Phase Extraction



Three Goals of SPE

1. Purification/Clean up

Before SPE




After SPE



Processing Samples

with a Vacuum Manifold, a pump and a waste trap

Sample reservoir is connected to cartridge by a tubing and adaptor cap. Atmospheric pressure pushes water sample through cartridge and flow is controlled by stopcocks between cartridge and cover of manifold.



SPE of Aqueous Samples

Example: A green sample containing a blue (non-polar) and yellow (polar) dye can be separated with SPE.



Sample Matrix

green dye



Contaminants

yellow = hydrophilic
(polar)



Target Analyze(s)

blue = hydrophobic
(non-polar)

Steps of SPE



1a. Condition

The two steps remove air and wet phase; phase is extended to its max. level.



1b. Equilibrate



2. Load sample

Analytes are retained by the sorbent.



3. Wash

Contaminants are washed off the phase.



4. Elute

Analytes are eluted from the phase.

2. Retention Mechanisms / Sorbent Selectivity



strata™ X

phenomenex®
...breaking with tradition™

SPE of Analytes from Water

Analytes in Water

SPE Phase must have a high retentivity towards analytes.
Interactions must be stronger than analyte's solubility in water.




Retention Mechanism

Hydrophobic / Non-Polar / Reversed Phase

Typical phases: C18, C8, Phenyl, Polymers (i.e. PSDVB, strata-X)

Most common for aqueous samples / Huge difference in polarity between phase and water.
Organic compounds interact with organic solid phase due to their organic frame.

Hydrophobic / Non-Polar / Reverse Phase
Interactions with C8 Functional Groups




Retention Mechanism

Hydrophilic / Polar / Normal Phase

Typical phases: Silica, Florisil[®], NH₂, CN

In most cases used for strong organic and non-aqueous environmental samples like hexane.

Hydrophilic / Polar / Normal Phase
Interactions with Si and CN
Functional Groups



↔ Hydrogen bonding and Dipole - Dipole interactions


Retention Mechanism

Ion Exchange – Ionic interactions

Typical phases: SCX, WCX, SAX, WAX (NH₂), strata-X-C, -X-CW, -X-AW

Can be used for aqueous environmental samples if just a specific ionic class of analytes has to be extracted. Adjusting appropriate pH for load and elution step is necessary.

Ionic / Ion Exchange Interactions with SCX and WCX Functional Groups



Matrix

Polarity of sample and solid phase must be opposed.

Liquid Phase / Sample	Solid Phase
Aqueous: Water, biological fluids, aqueous homogenates	Non-Polar C18, C8, Phenyl, SDB, strata-X
Organic: Non-polar organic solvents	Polar Silica, Florisil®, NH ₂ , CN
Aqueous or Organic	Ion Exchange SCX, WCX, SAX, WAX, strata-X-C, -X-CW, -X-AW

Sorbent Selectivity

Selection of sorbent dependent on analyte and sample matrix.
These sorbents have a single retention mechanism.

Most common
RP phases { **SDB Polymer**

C18E

C18U

C8

**Classic
Sorbents**

Phenyl

CN


NH₂

Most common
NP phases {

Si

FL-PR

SCX, WCX, SAX



Non-polar

Reversed Phase Sorbents
for organic compounds in
aqueous environmental
samples, e.g. surface water.

Polar

Normal Phase Sorbents for
polar organic compounds in
organic solvents e.g. organic
soil extracts.

Ionic

IEX Sorbents for specific
ionic compounds from
aqueous or organic samples.

Sorbent Selectivity

Question: Can a sorbent have more than one retention mechanism?

Multiple Modes

Mixed Phases: Hydrophobic + Ion exchange

Typical phases: strata Screen C and strata Screen A

“Universal”: Hydrophobic + Hydrophilic

Typical phases: modified polymers such as strata-X

“Universal”: Hydrophobic + Hydrophilic + Ion exchange

Typical phases: modified polymers such as strata-X-C(W), -AW

strata™ X

phenomenex®
...breaking with tradition™

3. strata-X for Water Analysis

strata[™]X

phenomenex[®]
...breaking with tradition[™]

Modern Polymer Sorbents

strata-X / strata-X-C(W) / strata-X-AW

strata-X: hydrophobic and hydrophilic retention mechanism
for neutral non-polar and polar compounds

strata-X-C: strata-X + strong cation exchange for bases

strata-X-CW: strata-X + weak cation exchange for strong bases

strata-X-AW: strata-X + weak anion exchange for acids

33 μm and 100 μm particles with large surface area $> 800\text{m}^2/\text{g}$ / $500\text{m}^2/\text{g}$

pH stability 1-14

strata™ X

phenomenex®
...breaking with tradition™


Advantages of strata-X Series

- Simplicity: One method approach
- High & reproducible recoveries for acids, bases and neutrals
⇒ Can be used for all kinds of environmental aqueous samples
- High pH stability (pH 1-14)
- High capacity → lower solvent consumption
- Resistant to de-conditioning
- Clean = LC/MS and GC/MS friendly

strata™ X

phenomenex®
...breaking with tradition™

strata-X simplifies Method Development



Step 1

Condition 5mL
of methanol.
Equilibrate with
5mL water

Step 2

Load sample

Step 3

Wash impurities
with 5mL of a 5%
methanol solution
in water

Step 4

Elute analyte
with 3-5mL of
methanol!

Method based on a 200mg sorbent bed mass

strata™ X

phenomenex®
...breaking with tradition™

4. Applications of Water Analysis

strata[™]X

phenomenex[®]
...breaking with tradition[™]

Multiresidue solid phase extraction of 16 pesticides from groundwater coupled with high-performance liquid chromatography *

SPE Method

Sorbent: strata-X, 200mg/6mL

Conditioning: 10mL dichloromethane

10mL acetonitrile

10mL water

Load: 1L water sample

Wash: water/methanol (95:5 v/v)

Elution: 5mL acetonitrile

5mL methanol

Pesticide	R%	Strata X	
			LOQ ($\mu\text{g l}^{-1}$)
Desisopropylatrazine		82(2)	0.009
Desethylatrazine		92(3)	0.009
Aldicarb		71(4)	0.14
Simazine		85(5)	0.009
Carbofuran		97(5)	0.06
Metalaxyl	100(4)	0.06
Atrazine		90(4)	0.009
2,4-D		113(4)	0.07
Metazachlor		90(4)	0.06
Dicloran		84(5)	0.03
Phenmedipham		66(2)	0.05
Linuron		93(2)	0.05
Iprodione		71(3)	0.05
Procymidone		74(6)	0.07
Fenitrothion		97(1)	0.15
Vinclozolin		88(10)	0.1

* Talanta 71 (2007) 25–30, Angelo Antonio D'Archivio, Maria Fanelli, Pietro Mazzeo, Fabrizio Ruggieri
*Universit'a degli Studi di L'Aquila, Dipartimento di Chimica, Ingegneria Chimica e Materiali,
 Via Vetoio, 67010 Coppito, L'Aquila, Italy*

strata™ X

phenomenex®
 ...breaking with tradition™

Solid phase extraction of pesticides from surface water *

SPE Method

Sorbent: strata-X, 200mg/6mL

Conditioning: 3x 1mL MeOH/ACN (1:1, v/v)

3x 2mL water

Load: 1L water sample (+ 0.5g NaCl)

Dry: 1h nitrogen stream

Elution: 3x 1mL MeOH/ACN (1:1, v/v)

Pesticides	Recovery [%]	RSD [%] n = 4	Set value [ng/l]
Desethylsimazin	85	5,8	100
2,6-Dichlorbenzamid	84	5,9	100
Ethidimuron	86	4,9	200
Chloridazon	86	3,6	50
Desethylatrazin	87	6,5	100
Desethylsebuthylazin	84	4,9	100
Bromacil	77	5,1	200
Simazin	84	7,1	100
Metribuzin	76	7,9	200
Desethylterbuthylazin	82	6,7	100
Metabenzthiazuron	82	5,4	50
Chlortoluron	83	4,5	50
Atrazin	82	9,0	100
Diuron	81	5,6	100
Isoproturon	81	5,5	100
Metazachlor	78	6,7	100
Terbumeton	77	9,1	100
Sebuthylazin	83	5,5	50
Propazin	79	10,3	50
Dimefuron	79	1,9	100
Terbuthylazin	78	9,0	100
Triadimenol	84	0,9	100
Epoxiconazol	80	2,6	200
Terbutryn	77	6,3	100
Metolachlor	71	16,2	200
Propiconazol	81	4,1	200
Kresoxim-methyl	79	3,5	100

* By courtesy of Mr. Cornelius Schmidt-Leistner from the „Bayerisches Landesamt für Wasserwirtschaft“ in Munich.

strata™ X

phenomenex®
...breaking with tradition™

Determination of selected human pharmaceutical compounds in effluent and surface water samples by high-performance liquid chromatography *

SPE Method

Sorbent: strata-X, 200mg/6mL

Conditioning: 3x 2mL methanol

3x 2mL water

3x 2mL water, pH 3

Load: 1L water sample, pH 3

Dry: vacuum aspiration

Elution: 3x 2mL methanol

Compound	SPE column recoveries ^a	
	Strata X R.S.D. (%)	LOD ^b (ng l ⁻¹)
Acetyl-sulfamethoxazole	56 (5.4)	50
Clofibric acid	83 (7.0)	50
Dextropropoxyphene	63 (3.9)	20
Diclofenac	62 (20)	20
Erythromycin	73 (30)	10
Ibuprofen	117 (22)	20
Lofepamine	4.2 (35)	10
Mefenamic acid	24 (7.9)	50
Paracetamol	75 (6.9)	50
Propranolol	45 (5.6)	10
Sulfamethoxazole	120 (16)	50
Tamoxifen	42 (40)	10
Trimethoprim	123 (2.5)	10

* Journal of Chromatography A, 1015 (2003) 129–141, Martin J. Hilton, Kevin V. Thomas
Centre for Environment, Fisheries and Aquaculture Science (CEFAS), Burnham Laboratory,
Remembrance Avenue, Burnham on Crouch, Essex CM0 8HA, UK

strata™ X

phenomenex®
...breaking with tradition™

Quantification of β -lactam antibiotics in urine and wipe samples from environmental and biological monitoring by SPE and LC-MS/MS

Institute of Energy and Environmental Technology (IUTA), D-47229 Duisburg, Germany

Table 2 Retention times (RT), recoveries and limits of detection (LOD) for wipe sample extracts (water) and urine samples

	RT [min]	Recovery Strata X [%] ^a		LOD water [μ g/L]	LOD urine [μ g/L]
Cefotiamе (1)	2.02	75		27	5.0
Amoxicillin (2)	2.81	78		0.8	20
Ampicillin (3)	7.02	74		0.5	0.4
Cefazoline (4)	9.32	71		0.9	11
Cefuroxime (5)	9.87	80		0.3	2.2
Piperacillin (6)	11.79	63		0.2	2.8
Penicillin G (7)	12.26	81		0.2	0.7
Penicillin V (8)	13.02	67		0.3	0.4

^a n = 1; ^b n = 3

*Thank you
for your attention*



strata™X

phenomenex®
...breaking with tradition™

strata-X is a trademark of Phenomenex Inc. Florisil is a registered trademark of U.S. Silica Company.