

**Mallinckrodt  
Baker**



# Contemporary methods of solid phase extraction

*Ivana Ivanovic*

Product Manager Separation Technology

Novi Sad, Serbia, September 9th, 2009



***Mallinckrodt Baker***

## OUTLOOK

- **J.T.Baker Solid-Phase Extraction Products**
- **Freon replacement – ISO-9377-2 (2000)**
- **Environmental Testing – Application examples**



***Mallinckrodt Baker***

## Principles of spe

---

SPE is an extraction process whereby, an aqueous sample is filtered through a thin bed of sorbent particles, the analytes of interest are removed from the sample matrix, and concentrated onto the sorbent. Once concentrated, the analytes are removed by an eluting solvent.



***Mallinckrodt Baker***


## What are the Benefits of SPE?

- ✓ SPE uses less solvent than LLE
- ✓ SPE is faster (at least 5 times)
- ✓ High capacity (2-8 %)
- ✓ Total SPE costs are considerably less than LLE
- ✓ SPE typically provides more accurate data than LLE
- ✓ High selectivity: broad choice of bonded phases and solvents
- ✓ Automation much easier



**Mallinckrodt Baker**

## The 4 Basic Steps in SPE





***Mallinckrodt Baker***


## Solid-Phase Extraction Flash Presentation




**Mallinckrodt Baker**

## J.T. Baker SPE sorbents portfolio

“Dry packed” 40  $\mu\text{m}$



“Slurry packed”





**Mallinckrodt Baker**

# Sorbents and configurations for BAKERBOND spe

**spe<sup>TM</sup>**

**Speedisk<sup>®</sup>**

## • Silica



## • Polymeric








Pre-cleaned housing

**BAKERBOND *Speedisk*™**  
Products are Protected by  
U.S. Patent No. 5,595,653

***SPEEDISK*®**






**Mallinckrodt Baker**

## Benefits of *Speedisk*® configuration

### Laminar BAKERBOND *Speedisk*®

- High flow rate
- High capacity
- Minimum of clogging



40 µm  
Conventional  
flow pattern



10 µm  
*Speedisks*®  
flow pattern



**Mallinckrodt Baker**

## Sold-Phase Extraction

# Conventional speed VS *Speedisk*®



Sorbent Weight	100 mg (40 µm) / 1 cc	10 mg (10 µm) / 1 cc
Sample Volume	2.0 ml	1.0 ml
Conditioning solvent	2.0 ml (20 sec)	0.5 ml (5 sec)
Sample addition	2.0 ml (100 sec)	1.0 ml (50 sec)
Wash solvent	1.5 ml (15 sec)	0.4 ml (5 sec)
Eluting solvent	2.0 ml (20 sec)	0.3 ml (3 sec)
Evaporation/reconst.	3 – 10 min	not necessary




## ハイパフォーマンス固相抽出 — 従来のSPEとの比較 —

*Speedisk*®  
concept

### 1 Speedisk Column は高流速での処理が可能

従来品(SPE Column)では流速を上げるとバンド幅が大きく広がっているが SpeediskColumnにはその影響がほとんど見られない。  
(サンプル:メチルオレンジ 1mg/ml)




### 2 Speedisk Column の大きい保持量



メチルオレンジを共に5mlづつ添加した。SpeediskColumn はバンドがシャープな為、サンプル添加時にサンプルをロスする事なく多くのサンプルを保持できる。

### 3 Speedisk Column の高い濃縮効果

メチルオレンジを同量づつ添加し、メタノールで100μl刻みにメチルオレンジの色が落ちるまで溶出させた。







**Mallinckrodt Baker**

## Disk and Column type configuration



### *Speedisk® Disk type*

- ✓ Large Volume Samples
- ✓ Environmental Applications



### *Speedisk® Column type*

- ✓ Low Volume Samples
- ✓ Biological/Pharmaceutical Applications



**Mallinckrodt Baker**

## *Speedisk*® silica based product range

- C<sub>18</sub>, C<sub>18</sub> Light Load, C<sub>18</sub> Polar Plus
- C<sub>8</sub>, C<sub>8</sub> Polar Plus
- Phenyl
- C<sub>4</sub>
- C<sub>2</sub>
- Cyano (CN)
- Amino (NH<sub>2</sub>)
- Diol (COHCOH)
- CBx WP, PEI WP, Butyl WP, HI Propyl WP (biotechnology)
- Silica (SiOH)
- Quaternary Amine (N<sup>+</sup>)
- Aromatic Sulfonic Acid (C<sub>6</sub>H<sub>5</sub>-SO<sub>3</sub>H)
- Carboxylic Acid (COOH)
- narc-1, narc-2 (for drugs of abuse analysis)





**Mallinckrodt Baker**

## Some typical *Speedisk*® Application

### Environmental

- - PAH 's from water and soil
- - PCB's from oil
- - Explosives in soil
- - Pesticides from water/soil
- - Phenoxy acid herbicides from water

### Food/Feed/Beverages

- - Aflatoxine from corn meal
- - Caffeine from di-cafeinated diet cola
- - Vitamin E from juice

### Pharmaceutical/Clinical/Biological

- - Benzodiazepines from serum
- - Anabolic Steroids/Urine
- - Aflatoxine from liver




***Mallinckrodt Baker***




**Accessories....**

**BAKER 12 G PTFE DESIGN**



**BAKER 12 G REGULAR DESIGN**



**BAKER 24 G REGULAR DESIGN**





***Mallinckrodt Baker***



*Speedisk*® Multi-96 Well Plate

***Speedisk*® 48 Processor:**

- Positive Pressure Processor for processing 1, 3, 6 ml columns in batches of 1 - 48 samples.

***Speedisk*® 96 Processor:**

- Positive Pressure Processor for processing 1 ml rimless columns in a 96 (8 x 12) microplate format



**Mallinckrodt Baker**

## Special designed *Speedisk* extraction manifolds

### ➤ Expanded Extraction Station 6-port



# Selection Guide for SPE sorbents and solvents

## Organic samples mw < 2000 (in solution)

This Sorbent Selection Guide is a systematic guide classifying sorbents according to polarity, ionizability and solubility in water or organic solvents. This information is useful in selecting the necessary components of a preliminary extraction method.

SOLUBLE IN	Organic Solvent Soluble			Water Soluble				
	Organic Polar	Organic Moderately Polar	Aqueous Non Polar	Ionic		Non-Ionic/non-paired		
				Aqueous Anionic	Aqueous Cationic	Aqueous Non Polar	Aqueous Moderately Polar	Aqueous Polar
MECHANISM	NPC	LSC	RPC	IEC	EC	RPC	LSC	NPC
SPE PHASE	H <sub>2</sub> O-Phobic DVB Cyan Diol Amino 1,2-Amino	H <sub>2</sub> O-Phobic DVB H <sub>2</sub> O-Phobic DVB Silica gel Florisil Alumina	H <sub>2</sub> O-Phobic DVB H <sub>2</sub> O-Phobic DVB SCX (NDC-2) Octadecyl Cyclodextrin Phenyl Dyads	H <sub>2</sub> O-Phobic WA-DVB H <sub>2</sub> O-Phobic SA-DVB Amino 1,2-Amino Quaternary Amino	H <sub>2</sub> O-Phobic SC-DVB H <sub>2</sub> O-Phobic SC-DVB Cyan Carboxylic Acid Sulfonic Acid	H <sub>2</sub> O-Phobic DVB H <sub>2</sub> O-Phobic DVB SCB-1/SCB-2 Octadecyl Cyclodextrin Phenyl Dyads	H <sub>2</sub> O-Phobic DVB H <sub>2</sub> O-Phobic DVB Silica gel Florisil Alumina	H <sub>2</sub> O-Phobic DVB Cyan Diol Amino 1,2-Amino
SOLVENTS*	Hexane Chloroform Dichloromethane Acetone Methanol	Hexane Chloroform Dichloromethane Ethyl acetate Methanol	Hexane Dichloromethane Acetone Acetonitrile Methanol Water	Acids, buffers	Acids, bases, buffers	Hexane Dichloromethane Acetone Acetonitrile Methanol Water	Hexane Chloroform Dichloromethane Ethyl acetate Methanol	Hexane Chloroform Dichloromethane Acetone Methanol

### \*Separation Mechanisms

LSC: Liquid Solid Chromatography (Adsorption)

NPC: Normal Phase Chromatography (Bonded Phase Partition)

RPC: Reversed Phase Chromatography (Bonded Phase Partition)

IEC: Ion-Exchange Chromatography (Bonded Phase Ion-Exchange)

### \*Sorbent phases listed in order of increasing polarity

\*Solvent solvents listed in order of increasing polarity

\*Solvent elution can be performed by combining two or more suitable solvents in order to achieve various degrees of polarity

### \*\*Solvents

0262 Hexane, ULTRA RESI-ANALYZED\*

0267 Chloroform, ULTRA RESI-ANALYZED\*

0264 Dichloromethane, ULTRA RESI-ANALYZED\*

0265 Ethyl acetate, ULTRA RESI-ANALYZED\*

0254 Acetone, ULTRA RESI-ANALYZED\*

0255 Acetonitrile, ULTRA RESI-ANALYZED\*

0077 Methanol, ULTRA RESI-ANALYZED\*

42 9 Water, ULTRA RESI-ANALYZED\*

SCX: Styrene Divinyl Benzene

DVB: Divinyl Benzene

H<sub>2</sub>O-Phobic WA-DVB: Weak anion exchange

H<sub>2</sub>O-Phobic SC-DVB: Strong cation exchange

H<sub>2</sub>O-Phobic SA-DVB: Strong anion exchange

H<sub>2</sub>O-Phobic SC-DVB: Strong cation exchange

This Sorbent Selection Guide is also available on:

[www.jtbaker.com/chromatography/SolidPhaseExtraction.asp](http://www.jtbaker.com/chromatography/SolidPhaseExtraction.asp)

[www.jtbaker.com](http://www.jtbaker.com)

**tyco** / **Huntsman** / **Millipore** / **Baker**

Baker Instrument Sales  
Department 25  
P.O. Box 1  
10000 Cypress  
Houston, Texas 77036

Baker Instrument Sales  
c/o J.T. Baker Ltd.  
London  
E10 5JF  
United Kingdom

Baker Instrument Sales  
2000 Highway 100  
P.O. Box 100  
Bly, Georgia 30516

Baker Instrument Sales  
Innovative Systems  
10000 Cypress  
Houston, Texas 77036

Baker Instrument Sales  
2000 Highway 100  
P.O. Box 100  
Bly, Georgia 30516

Baker Instrument Sales  
2000 Highway 100  
P.O. Box 100  
Bly, Georgia 30516

**A PART OF**  
**ANALYTICAL**  
**LABORATORY**



**Sample Preparation**  
When every step counts

**Conventional SPE** **JTBaker Speedisk**



**Speedisk™ columns:**  
**9 times faster** than conventional SPE.

Now you can sample in just a few minutes. The sample preparation process is faster, easier, and more efficient. The JTBaker Speedisk™ column is the only one that can be used for both solid and liquid samples. The JTBaker Speedisk™ column is the only one that can be used for both solid and liquid samples. The JTBaker Speedisk™ column is the only one that can be used for both solid and liquid samples.



***Mallinckrodt Baker***



Development of new products  
applicable to  
the quantitative determination  
of the mineral oil index (hydrocarbon index)  
according to ISO-9377-2 (2000)





**Mallinckrodt Baker**

# Determination techniques for Mineral Oil in Europe



## Use of Infrared (IR) since 1951


NEN 6673 (1981), DIN 38409-H18 (1981), NEN 6675 (1989)

Tetrachloromethane > Freon > Tetrachloroethylene

## From 1997 Gaschromatographic determination

NEN 5733 (1997), DIN-EN ISO 9377-4 (1999), **ISO 9377-2 (2000)**

Petroleum Ether, Iso-Hexane, Pentane, n-Hexane etc.



**Table 2.** Comparison of chlorofluorocarbon replacements to trichlorofluoromethane.

Compound	Intended use	ODP	GWP	Estimated atmospheric lifetime, years
Trichlorofluoromethane	Refrigeration, blowing agent	1	1	144
1,1-Dichloro-1-fluoroethane (HCFC-141b)	Blowing agent, solvent	0.1	0.12	11
1,1,1,2-Tetrafluoroethane (HFC-134a)	Refrigeration	—	0.4	14
Pentafluoroethane (HFC-125)	Refrigeration, fire extinguisher	—	0.84	41
1-Chloro-1,2,2,2-tetra-fluoroethane (HCFC-124)	Blowing agent, refrigeration	0.02	0.1	1.3
1,1,1-Trifluoro-2,2-dichloroethane (HCFC-123)	Refrigeration, blowing agent	0.02	0.02	1.6

Abbreviations: ODP, ozone depletion potential; GWP, global warming potential.




**Mallinckrodt Baker**

## EXTRACTION SOLVENTS (IR-determination)

### FREON


(1,1,2-Trichlorotrifluoroethane)



Freon is used as an extraction agent for mineral oil determination using IR-detection since 1951. Although new methods are developed utilizing a sophisticated GC apparatus, still a lot of laboratories deal with a straightforward IR analyser. Remark: Freon is almost fully prohibited from the European Market

### Tetrachlorethylene

(Perchloroethylene) *product code 9360*



Tetrachlorethylene is a largely used alternative for Freon. This solvent also acts as an extraction agent for the determination of the Hydrocarbon oil index, using IR detection. J.T. Baker's Ultra Resi Tetrachloroethylene suits to your existing Freon application.



**Mallinckrodt Baker**

## **ISO-9377-2 (2000)**

**Since October 2000 an alternative method, to replace Freon,  
has been directed by the EC (ISO-9377-2)**

- Single hydrocarbon or a mixture of hydrocarbons, both with a Boiling range between 36-69°C, replaces Freon
- The extraction is followed by cleanup for removal of polar substances (using Activated Florisil®) and Gas Chromatographic (GC) separation with on-column injection (including large volume injection (LVI) technique) and Flame Ionization Detection (FID).





**Mallinckrodt Baker**



## ISO-9377-2

### Hydrocarbon oil index

Sum of concentrations of compounds extractable with a hydrocarbon solvent with boiling point between 36 °C and 69 °C, not adsorbed on Florisil® and which may be chromatographed with retention times between those of n-Decane (C<sub>10</sub>H<sub>22</sub>) and n-Tetracontane (C<sub>40</sub>H<sub>82</sub>)



**Mallinckrodt Baker**


## Procedure ISO 9377-2 (2000)- Summary



- **Water sample (1L) is extracted with 50 ml extracting agent**
  - ‘Single Hydrocarbon Solvent with boiling range 36 °C - 69 °C’
- **Water is removed with Sodium sulfate, anhydrous**
  - ‘ $\text{Na}_2\text{SO}_4$  has to dried before use’
- **Polar substances are removed by clean-up on Florisil**
  - ‘Florisil has to be activated by heating to 140 °C for 16 h’
- **Preconcentration step with evaporation apparatus**
  - ‘To enhance the limit of detection’
- **Final analysis of hydrocarbon index with on-column GC-FID / Large Volume Injection Technique**



**Mallinckrodt Baker**



Ultra Trace GC Large Volume Injection of Iso-Hexane inclusive C10/C40 internal standards, injection volume 100  $\mu$ l (by kind permission of Interscience B.V.)



**Mallinckrodt Baker**

## Products for Mineral Oil Determination according to ISO 9377-2 (2000)



9266 Petroleum Ether

9267 Iso-Hexane

- *total peaks between n-decane ( $C_{10}H_{22}$ ) and n-tetracontane ( $C_{40}H_{82}$ )...max.  
0.5 mg/L*

7061-00 Florisil, activated

- *(16 h at 140 °C)*

3377-00 Sodium Sulfate, anhydrous

7495-18 SPE Column Ready to use with 2 g Sodium Sulfate / 2 g Florisil

7495-04 SPE Column Ready to use with 0.5 g Sodium Sulfate / 0.5 g Florisil

*For use in Large Volume injection techniques*

0252 Sand "washed and ignited"

*Used as blank sample*

0168 Magnesium Sulfate Heptahydrate

*Agent to avoid formation of emulsions*

Info Bulletin



Adobe Acrobat  
7.0 Document



***Mallinckrodt Baker***

## Application Examples – Environmental

- ◆ Pesticides (Triazines, Organchlorines etc)
- ◆ Polychlorinated Biphenyls (PCBs)
- ◆ Polycyclic Aromatic Hydrocarbons (PAH's)
- ◆ Acrylamide
- ◆ Pharmaceuticals
- ◆ Endocrine disruptors



**Mallinckrodt Baker**

# SPE Application Notes Examples-Environmental

**Extraction and determination of pharmaceutical residues and related polar contaminants in water samples**



Adobe Acrobat  
7.0 Document

**Clean up procedure with Florisil®/Sodium Sulfate for Hydrocarbon Index determination according to ISO 9377-2.**



Adobe Acrobat  
7.0 Document

**Multiresidue Analysis Method of Triazines, Organochlorine Pesticides and Polyaromatic Hydrocarbons in Drinking Water (adapted from the method E.P.A. 525.2.)**



Adobe Acrobat  
7.0 Document

**Extraction of Polycyclic Aromatic Hydrocarbons in water**



Adobe Acrobat  
7.0 Document

**Extraction of steroid hormones, hormone conjugates and macrolide antibiotics in influents and effluents of sewage treatment plants**



Adobe Acrobat  
7.0 Document



***Mallinckrodt Baker***

## SPE Application Lists



Microsoft Word  
Document

**SPE and Speedisk Applications**



Microsoft Word  
Document

**Deventer Applications**



Microsoft Word  
Document

**Evaluation Purpose Only**



Microsoft Word  
Document

**Zanob Applications**





***Mallinckrodt Baker***

# QUESTIONS?







***Mallinckrodt Baker***

**THANK YOU FOR YOUR  
ATTENTION!**





***Mallinckrodt Baker***

Solid-phase extraction

**BAKERBOND™**

**Activated Carbon**





**Mallinckrodt Baker**

**Solid-phase extraction**

## **BAKERBOND™ ACTIVATED CARBON PRODUCTS**


- A. BAKERBOND™ Activated Carbon**
- B. specially developed for the extraction of Acrylamide from surface and drinking water**
- C. German standard method (DIN 38 413 - P xx)\* for the determination of acrylamide in water (e.g. drinking water, surface water) by high performance liquid chromatography with mass spectrometric detection (HPLC MS/MS)**



**Mallinckrodt Baker**

Solid-phase extraction

## Acrylamide



- ▶ One of the contaminants in 1<sup>st</sup> Priority List of FDA, WHO and EPA
- ▶ Formed during the cooking of carbohydrate rich food



***Mallinckrodt Baker***

## Solid-phase extraction

- ▶ Fried, baked, or roasted food
- ▶ Present in breads and cereals



- ▶ Preparation of polyacrylamide or its copolymers
- ▶ Polyacrylamide used in purification of drinking, waste and process water





**Mallinckrodt Baker**

**Solid-phase extraction**

## PHYSICAL PROPERTIES OF BAKERBOND™ ACTIVATED CARBON

<b>BAKERBOND™ Carbon</b>	
<b>Appearance</b>	<b>Black Spheres</b>
<b>Particle Size</b> (by Laser Scanning in Water Suspension) d <sub>10</sub> [μm] d <sub>50</sub> [μm] d <sub>90</sub> [μm]	210 – 240 300 – 330 405 – 435
<b>Specific Surface Area</b> (by Nitrogen Adsorption (BET-method), m <sup>2</sup> /g)	<b>Typical 1300 !!!</b>
<b>Porosity</b> <b>Microporosity (cm<sup>3</sup>/g)</b> <b>Mesoporosity (cm<sup>3</sup>/g)</b>	typical 0.75 typical 0.50
<b>Loss on drying (%)</b>	1.40 – 1.60
<b>Bulk density (g/ml)</b>	typical 0.40



**Mallinckrodt Baker**

Solid-phase extraction

## BAKERBOND™ Activated Carbon

### Sorbent

PN 7532-00 – 100 gr of Activated Carbon



Microsoft Word  
Document

Sales Tools

### SPE columns

PN 7575-06 – 6 ml PP / 500 mg

PN 7575-07 – 6ml PP / 1000 mg



Microsoft Word  
Document

Organophosphorous Pesticides

Application



Microsoft Word  
Document

Acrylamide Application



Microsoft Word  
Document


Acidic Herbicides Application



**Mallinckrodt Baker**

Solid-phase extraction

## BAKERBOND™ Carbon-Amino Double Phase Column



**Clean-up of Pesticides and Insecticides from Food and Feed**



Microsoft Word  
Document

**Sales Tool**



Microsoft Word  
Document

**Application**

**Pesticide residue analysis in agricultural products**  
**Organophosphates, Carbamates, Pyrethroids**





**Mallinckrodt Baker**

**Solid-phase extraction**

## **BAKERBOND™ SDB-1**

**Extraction of Chlorophenoxyacid Herbicides from Water**

**Extraction of Phenols from Water**

**Extraction of Carbamates from Water**

**Extraction of Explosives and Metabolites from Water**

**Extraction of Organophosphorus Pesticides from Water**

**Extraction of Pesticides from Water using**

**Extraction of Pesticides from Soil and Urine**



Microsoft Word  
Document

**Sales Tool**



Microsoft Word  
Document

**Instructions for  
use**



Microsoft Word  
Document

**ISO norm : *Final Draft ISO/CD 18857-2 „Water quality – Determination of selected alkylphenols, alkylphenol ethoxylates and bisphenol A-Method for non filtered samples using solid-phase extraction and gas chromatography with mass selective detection - to be ready soon!***



**Mallinckrodt Baker**

## SPE Application Examples – Pesticides

Pesticide class	ED Chemicals in the Group
<b>Carbamates</b>	Aldicarb, Carbaryl, Methomyl, Baygon (propoxur), Bendiocarb, Oxamyl
<b>Organochlorines</b>	Aldrin, Chlordane, Endosulfan, Endrin, Nonachlor, Oxychlordane
<b>Linuron, diuron, and derivatives/metabolites</b>	Ethylene Thiourea (ETU), Linuron (Lorox) Diflubenzuron, Diuron
<b>Organophosphates</b>	Malathion, Ethylmalathion, Methylmalathion, Chlorpyrifos, Acephate, Chlorfenvinphos, Diazinon, Dichlorvos, Dimethoate, Fenthion, Glufosinate, Mevinphos, Parathion, Phosphamidon, Quinalphos,
<b>Triazines and triazoles</b>	Amitrol, Atrazine, Biteranol, Cyanazine, Simazine, Terbutryn, Triadimefon, Triadimenol, Triazines



**Mallinckrodt Baker**

## SPE Application Examples – Pesticides



HTML Document

**EXTRACTION OF EPA METHOD 8081A OR 8082 ANALYTES  
ORGANOCHLORINE PESTICIDES OR POLYCHLORINATED BIPHENYLS  
USING H<sub>2</sub>O-PHOBIC DVB EXTRACTION DISK (GC-ECD)**



HTML Document

**EXTRACTION OF ORGANOCHLORINE PESTICIDES, HERBICIDES, AND  
ORGANOHALIDES FROM DRINKING WATER EPA METHOD 508.1 (*Speedisk*  
C<sub>18</sub> High Capacity Extraction Disk, GC-ECD)**



HTML Document

**EXTRACTION OF PESTICIDES FROM WATER USING SDB-1 or  
*Speedisk*® H<sub>2</sub>O-Phobic DVB Column (HPLC)**




HTML Document

**EXTRACTION OF PESTICIDES FROM WATER USING C<sub>18</sub> POLAR PLUS®  
(GC or HPLC)**




**Mallinckrodt Baker**


# SPE Application Examples – Pesticides




**EXTRACTION OF CARBAMATES FROM WATER USING SDB-1  
OR *Speedisk*® H<sub>2</sub>O-Phobic DVB Column (HPLC)**



**EXTRACTION OF TRIAZINE HERBICIDES FROM SOIL  
(SIMAZINE, ATRAZINE, PROPAZINE) – (Aromatic Sulfonic acid  
SPE, HPLC)**



**EXTRACTION OF TRIAZINES AND URONES FROM WATER  
USING C18 POLAR PLUS® (HPLC)**








**EXTRACTION OF TRIAZINES FROM WATER USING C18  
POLAR PLUS® (HPLC)**



**Mallinckrodt Baker**

# PolyChlorinated Biphenyl's – PCBs

- Class of organic compounds characterized by two benzene rings linked by a C-C bond
- Very resistant to degradation (oxidation, acids, bases, temperature)



- Soluble in most of the common organic solvents, slightly soluble in water
- Used as cooling/isolation fluids and fire retardants
- Accumulated - high stability under environmental conditions



**Mallinckrodt Baker**

# SPE Application Notes Examples – PCBs



HTML Document

**EXTRACTION OF POLYCHLORINATED BIPHENYLS (PCBs) FROM WATER– *Speedisk C<sub>18</sub>* (GC/MS determination)**



HTML Document

**EXTRACTION OF POLYCHLORINATED BIPHENYLS (PCBs) FROM WASTEWATER EPA Method 1668– *Speedisk C<sub>18</sub>* (GC/MS determination)**




HTML Document

**EXTRACTION OF PCBs FROM TRANSFORMER OIL – Florisil (GC/MS determination)**





# Polycyclic Aromatic Hydrocarbon's – PAHs



- (1) Phenanthrene
- (2) Anthracene
- (3) Pyrene
- (4) Benz[a]anthracene
- (5) Chrysene
- (6) Naphthacene
- (7) Benzo[c]phenanthrene
- (8) Benzo[ghi]fluoranthene
- (9) Dibenzo[c,g]phenanthrene
- (10) Benzo[ghi]perylene
- (11) Triphenylene
- (12) o-Tephenyl
- (13) Benzo[a]pyrene
- (14) p-Tephenyl
- (15) Benzo[a]pyrene
- (16) Tetrabenzonaphthalene
- (17) Phenanthro[3,4-c]phenanthrene
- (18) Coronene



**Mallinckrodt Baker**

# Polycyclic Aromatic Hydrocarbon's – PAHs

- ↳ Large and heterogeneous group of organic contaminants
- ↳ Formed and emitted as a result of the combustion of organic material
- ↳ Lipophilic compounds ( high affinity for organic matter)
- ↳ Differ substantially in their physicochemical properties
- ↳ Physico-chemical properties largely determine the environmental behavior



**Mallinckrodt Baker**

## SPE Application Notes Examples – PAHs



HTML Document

**EXTRACTION OF POLYCYCLIC AROMATIC HYDROCARBONS (PAHs)  
FROM WATER – *Speedisk C18* (GC/MS determination)**



HTML Document

**EXTRACTION OF POLYCYCLIC AROMATIC HYDROCARBONS (PAHs)  
FROM WATER - *Speedisk C<sub>18</sub>* High Capacity (GC/MS determination)**



HTML Document

**EXTRACTION OF POLYCYCLIC AROMATIC HYDROCARBONS (PAHs)  
FROM DRINKING WATER EPA Method 550.1 – *Speedisk C<sub>18</sub>* High Capacity  
(HPLC-UV/FLU determination)**



**Mallinckrodt Baker**

## Application Examples – Environmental

SORBENT	APPLICATION
C18	PAH's, Phthalates, Organochlorine pesticides, PCB's etc.
C18 XF	Extra filter for dirty samples
C18 Polar Plus	Phenols, Chlorophenoxy acids, Urones etc.
C8	Diquat, Paraquat
SAX	Haloacetic acids, Dalapon
H <sub>2</sub> O Phobic-DVB (HC)	Chlorinated Acids,
H <sub>2</sub> O Philic-DVB	Carbamates, Pharmaceutical residues from water
Oil & Grease	EPA 1664, Rev.A


XF= Extra filter

HC= High Capacity



**Mallinckrodt Baker**

## J.T.Baker Analytical BU program





***Mallinckrodt Baker***

# QUESTIONS?

